



KLINGER  
EDUCATIONAL  
PRODUCTS CORP.

112-19 14TH ROAD  
COLLEGE POINT, NEW YORK 11356  
(718) 461-1822

### The Hall effect for silver

Verifying the proportional relationship between Hall voltage  $U_H$  and magnetic flux density  $B$ , determining the polarity of the charge carriers which are mainly responsible for charge transport in silver. Calculating the charge carrier density  $n$ .

If a current-carrying metallic conductor band is in a magnetic field perpendicular to the direction of current, a transverse electrical field and a potential difference are produced (Hall effect).

The following equation is valid for the Hall voltage  $U_H$ :

$$U_H = \frac{1}{n \cdot e} \cdot \frac{B \cdot I}{d} \quad (1)$$

$B$ : Magnetic flux density

$I$ : Current

$d$ : Thickness of the band-shaped conductor

$e$ : Elementary charge

$n$ : Concentration of charge carriers

$\frac{1}{n \cdot e}$  is called the Hall constant  $R_H$ .  $R_H$  is a value which is dependent on material and temperature:

$$R_H = \frac{1}{n \cdot e} \quad (2)$$

The conductor band consists of silver in the following experiment. First, it will be proved that  $U_H \sim B$ .

The polarity of the charge carriers which are mainly responsible for the current (polarity of the Hall constants) can be determined from the direction of the Hall voltage. The concentration of the charge carriers is established experimentally since all the values in "1" except for  $n$  can be measured. The Hall voltage  $U_H$  is caused by deflection of the moving charge carriers in the magnetic field by the Lorentz force, whose direction may be predicted by the "right-hand rule".

### Apparatus:

1	Hall effect apparatus, silver .....	586 81
1	U-core with yoke .....	562 11
1	Clamping device .....	562 12
1	Pair of bored pole pieces .....	560 31
2	Coils, 250 turns .....	562 13
1	Extra-low variable transformer .....	522 39
2	Interchangeable scale demonstration meters, basic units .....	530 50
1	Measuring module 0 to 1000 mT/3000 mT .....	530 75
1	Measuring module 0-3 A/10 A .....	530 65
1	Shunt resistor 30 A .....	530 90
3	Connecting leads, red, 50 cm .....	501 25
2	Connecting leads, blue, 50 cm .....	501 26
3	Connecting leads, black, 50 cm .....	501 28
1	Connecting lead, black, 25 cm .....	501 23
1	Microvoltmeter .....	532 13
1	Pole probe .....	516 501
1	Power supply unit 0-12 V/20 A, regulated .....	522 47

### Recommended:

1	Calibrating magnet for pole probe ...	516 53
1	Multimeter 2H (15 A -) .....	531 53
2	Batteries 1.5 V/IEC R6 .....	685 44




Fig. 1: Experiment setup for the Hall effect.  
 $I_B$ : Field current

$I_Q$ : Transverse current

### Setting up:

#### Note:

Only switch on circuits for a short time for a transverse current above 15 A or magnetic currents above 5 A because the connecting leads will heat up and the coils designed for 5 A will overload.

In the transverse circuit, use leads which permit a load of 20 A (e.g. connecting leads 501 20-29 or safety connecting leads 500 65-74).

Set up the apparatus as in Fig. 1, initially without the Hall effect apparatus. Set up the electromagnet with exactly the pole piece spacing that is given by the thickness of the support plate of the Hall effect apparatus. To do this, loosen the clamping device and place one edge of the Hall effect apparatus between the pole pieces. Then push the latter as close as possible to the support plate.

Calibrate the Hall probe with the calibrating magnet as described in the instructions for use 516 501. Remove the protective sleeving from the pole probe.

### Carrying out the experiment:

#### a) $I_B$ - B calibrating curve:

Demagnetize the iron of the electromagnets before recording the I-B calibrating curve and before determining B from this curve (allow 5 A  $\sim$  to flow through the coils for a short time).


Measure the magnetic flux density B as a function of magnetic current  $I_B$ ; to do this, increase  $I_B$  in steps of 0.5 A (Fig. 3).

b)  $U_H$  is a function of B for constant transverse current  $I_Q$ :

Mount the Hall effect apparatus in the electromagnet as shown in Fig. 1, with the pole pieces pushed as close as possible to the support plate (air gap as narrow as possible and of the same width as for the recording of the calibrating curve).

Connect the Hall effect apparatus to the microvoltmeter and the power source as shown in Fig. 2. The field direction should be as printed on the support plate (current direction  $I_B$ ; coil connections).

Fig. 2: Electrical connections of the pole effect apparatus;  $I_Q$ : transverse current.



Reset the measuring instrument display for the Hall voltage  $U_H$  to zero before switching on the magnetic current  $I_F$ , but with the transverse current switched on. Correct the zero point by means of button (a) (Fig. 1) on the microvoltmeter.

After switching off the magnetic current, check the zero point again and, if necessary, take into account any fluctuations which have occurred.

Read off the respective zero voltage  $U_H$  for each

set current  $I_B$ . Read off the effective field strength from the  $I_B$ -B calibrating curve for each current  $I_B$ . Transverse current  $I_Q = 15$  A and  $I_Q = 20$  A (Fig. 4).

Polarity of the charge carriers:

To determine the polarity of the Hall voltage for the chosen current direction.

Determination of the charge carrier concentration  $n$  and Hall constant  $R_H$ :

Set a transverse current  $I_Q = 15$  A, and a field

current  $I_B = 8.5$  A ( $\approx B = 0.805$  Tesla, in

accordance with calibrating curve). Measure  $U_H$ .

Repeat measurement for  $I_Q = 20$  A,  $I_B = 8.5$  A.




Fig. 3:  $I_B$ -B calibrating curve; pole piece

distance the same as the thickness of the support plate

$I_B$ : Coil current

B: Magnetic flux density of the field produced by  $I_B$ . Saturation at large

field currents.

#### Measurement example:

a) Table 1/Fig. 3

Table 1:

I in A	B in T
0	0
0,5	0,118
1	0,200
1,5	0,295
2	0,374
2,5	0,455
3	0,520
3,5	0,585
4	0,630
4,5	0,665
5	0,695
5,5	0,715
6	0,735
6,5	0,748
7	0,760
7,5	0,780
8	0,790
8,5	0,800
9	0,810

b) Fig. 4




Fig. 4:  $U_H$  as a function of B; parameter  $I_Q$

B in accordance with calibrating curve in Fig. 3 from  $I_B$ ;

Hall foil made of silver.

c) The microvoltmeter indicates negative voltage values for the setup chosen in Fig. 2. The upper side of the Hall foil is thus charged negatively with respect to the lowest side.

d<sub>1</sub>):

$I_Q = 15 \text{ A}$ ;  $B = 0.805 \text{ T}$ ; foil thickness  $d = 5 \cdot 10^{-5} \text{ m}$ ;  
 $10^{-5} \text{ m}$ ;

$e = 1.602 \cdot 10^{-19} \text{ C}$ ;  $U_H = 1.87 \cdot 10^{-5} \text{ V}$ .

d<sub>2</sub>):

$I_Q = 20 \text{ A}$ ;  $B = 0.805 \text{ T}$ ;  $U_H = 2.4 \cdot 10^{-5} \text{ V}$ .

#### Evaluation and results:

For b)

The graphs in Fig. 4 show that  $U_H \sim B$  and that  $U_H$  increases with increasing transverse current  $I_Q$ .

Note:

Experimental proof of the proportional relationship  $U_H \sim I$  may easily be obtained by

measuring  $U_H$  for various  $I_Q$  (for constant field current  $I_B$ ).

For c)

Negative Hall voltages are obtained for the setup of current and field direction chosen in Fig. 2. If the "right-hand rule" is used, we find that the conduction mechanism for silver is mainly effected by negative charge carriers.

Note:

In 1916, Holman obtained certain proof that electrons are the charge carriers in metals.

For d<sub>1</sub>)

Evaluation in accordance with (1) and (2):

$$R_H = 7,74 \cdot 10^{-11} \text{ m}^3 \cdot \text{C}^{-1}$$

$$n = 8,06 \cdot 10^{28} \text{ m}^{-3}$$

For d<sub>2</sub>)

$$R_H = 7,45 \cdot 10^{-11} \text{ m}^3 \cdot \text{C}^{-1}$$

$$n = 8,37 \cdot 10^{28} \text{ m}^{-3}$$

Note:

$$\text{Theoretical value: } R_H = 8.9 \cdot 10^{-11} \text{ m}^3 \cdot \text{C}^{-1}$$

$$n = 6.6 \cdot 10^{28} \text{ m}^{-3}$$

$$(\text{atoms density } 5.8 \cdot 10^{22} \text{ cm}^{-3})$$